


COMPLETE DIAGNOSIS OF DRY EYE


4 Non-contact exams
meeting the criteria of DEWS 2*

LACRYDIAG™


INTERFEROMETRY




N.I.B.U.T.




TEAR MENISCUS



MEIBOGRAPHY




*TFOS DEWS II Diagnostic Methodology report - The Ocular Surface -Volume 15, Issue 3, July 2017, Pages 539-574
The full report is available online from: www.tfosdewreport.org




MEIBOGRAPHY

EVOLUTION OF MEIBOMIAN GLAND DYSFUNCTION



LACRYDIAG™

CLASSIFICATION INTO 5 LEVELS OF THE MEIBOMIAN GLAND DYSFUNCTION



For more information:
 Phone. +33 (0)4 73 745 745 - contact@quantel-medical.fr